

NEW!

SIMPLY SMARTER CRYOSURGERY



INTRODUCING



**CRYOSURGERY WITHOUT
THE GUESSWORK**

- ❄️ **Nitrous Oxide OR Carbon Dioxide Gas**—First and only cryo product that offers the flexibility to use either cryogen based on your patients' needs.
- ❄️ **Advanced Treatment Wand**—Now you can treat you patients in any direction, at any angle, on any part of the body.
- ❄️ **Preset Dispense Times**—Choose between 5-second or 10-second treatment times and minimize the risk of under- or over-treating.
- ❄️ **Gas Level Indicator**—Our new lighted indicator ensures you will know when to order a new tank.




CRYOSURGERY JUST GOT A WHOLE LOT SMARTER!

For the first time ever, you can choose between two cryogen—Nitrous Oxide (-89°C) or Carbon Dioxide (-79°C)—to get great first-time outcomes and minimize damage to surrounding skin.

N₂O is an ideal cryogen for most common lesions such as verruca vulgaris, angiomas, seborrheic keratosis and more.*

What cryogen is best for cosmetic lesions such as skin tags, age and sun spots on the hands and face? The answer is CO₂ which provides a more superficial treatment and minimizes your patients' downtime and recovery.



Item No.	Product Code	CryoLab® Medical Product Description
200-1000	CryoLab Medical Device Unit	<p>Kit Contains:</p> <ul style="list-style-type: none"> 1 Canister of cryogen 1 N₂O OR CO₂ 1 12V/20W Power Supply 1 90" Power Cord 1 Quick Start Guide 1 Instructions for Use One-year Limited Warranty*
000-2730	Replacement N2O Cylinder	<p>Contains:</p> <ul style="list-style-type: none"> 1 Canister of N2O cryogen 1 Bleed Plug 1 Return Instructions 1 Return Label
000-2731	Replacement CO2 Cylinder	<p>Contains:</p> <ul style="list-style-type: none"> 1 Canister of CO2 cryogen 1 Bleed Plug 1 Return Instructions 1 Return Label



After Action Medical and Dental Supply, LLC
 4444 Decatur Blvd Suite 100 • Indianapolis, IN 46241
afteractionmedical.com
 1-800-892-5352 • sales@afteractionmedical.com

*Please contact CryoConcepts, LP for details

*For complete list of indications, please see the Instructions for Use

